

The shaker is widely used in Western experiments such as blocking of blotting membranes, antibody incubation, washing, fixation, staining, and destaining of electrophoretic gels; it can also be used for shaking liquids in various centrifuge tubes, multi-well plates, Petri dishes, and culture bottles, mix well.

## Product Features:

1. Small design, easy to move
2. The pendulum mode can automatically stop to prevent the dry film from affecting the experimental results
3. Power-off memory, automatically continues to run after power-off



**MDS-100P**  
**Pendulum shake**



**MDS-100H**  
**Horizontal shake**

Digital Shaker	MDS-100P	MDS-100H
Movement	Pendulum Movement	Horizontal Movement
Speed range	0-200 rpm	0-200 rpm
Max load	2.5 Kg	2.5 Kg
Tray size	317*245mm	317*245mm
Dimension	317*245*139mm	317*245*130mm
Tilt Angle	20°	20°